

Product Datasheet

Radox - Interleukin-8 (IL-8) (RCP9125)

Description	Interleukin-8 (IL-8) belongs to the neutrophil-specific CXC family of chemokines. It is one of the i...
Preservatives	20 mM Tris-HCl pH 8.0, 250 mM NaCl, 50 % Glycerol
Form/Appearance	Liquid
Storage	at -20 °C
Tag	N terminal containing a HIS tag
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	13.7 kDa
Uniprot ID	P10145
Expiration Date	6 months from date of receipt.

