

Product Datasheet

Randox - C-Reactive Protein (CRP) (RCP9131)

Description	Human C reactive protein (CRP) is a member of the pentraxin family of serum proteins and is composed...
Preservatives	20 mM Tris-HCl pH 7.5, 0.14 M NaCl, 2 mM CaCl ₂ , 0.05% NaN ₃
Form/Appearance	Liquid
Storage	at 4 °C
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	27.18 kDa
Uniprot ID	P02741
Expiration Date	6 months from date of receipt.

