

Product Datasheet

Radox - Eotaxin (RCP9149)

Description	Eotaxin / CCL11 is a member of the beta (C-C) motif chemokine family of inflammatory and immunoregul...
Preservatives	25 mM Tris-HCl pH 8.0, 450 mM NaCl, 50 % glycerol
Form/Appearance	Liquid
Storage	at -20 °C
Tag	N terminal containing a HIS tag
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	12.86 kDa
Uniprot ID	P51671
Expiration Date	6 months from date of receipt.

