

Product Datasheet

Radox - C-Reactive Protein (CRP) (Western Blot Control) (RCP9293)

Description	Human C reactive protein (CRP) is a member of the pentraxin family of serum proteins and is composed...
Preservatives	1X Laemmli buffer
Form/Appearance	Liquid. Supplied in 1x Laemmli Buffer (25mM Tris-HCl pH 6.8, 50mM DTT, 1% (w/v) SDS, 0.1% (w/v) Bromophenol Blue, 2.5% Glycerol).
Storage	at 4 °C
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	27.18 kDa
Uniprot ID	P02741
Expiration Date	6 months from date of receipt.

