

## Product Datasheet

### Randox - C-Reactive Protein (CRP) (RCP9295)

<b>Description</b>	Human C reactive protein (CRP) is a member of the pentraxin family of serum proteins and is composed...
<b>Preservatives</b>	20 mM Tris-HCl pH 7.5, 0.14 M NaCl, 2 mM CaCl <sub>2</sub> , 0.05% NaN <sub>3</sub>
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	at 4 °C
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA, Western Blotting, etc.
<b>Purity</b>	> 95 %, by SDS-PAGE under reducing conditions
<b>Source</b>	E.coli -derived human recombinant protein
<b>MW</b>	27.18 kDa
<b>Uniprot ID</b>	<b>P02741</b>
<b>Expiration Date</b>	6 months from date of receipt.

