

Product Datasheet

Radox - Transforming Growth Factor Beta - 1 (TGF-B1) (RCP9395)

Description	TGF beta 1 is a member of the larger TGF beta superfamily of cytokines produced by a wide variety of...
Preservatives	25mM Sodium Acetate pH 4.8, 50% Glycerol
Form/Appearance	Liquid
Storage	at -20 °C
Tag	N terminal containing a HIS tag
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	17.3 kDa
Uniprot ID	P01137
Expiration Date	6 months from date of receipt.

