

Product Datasheet

Radox - Interleukin-33 (IL-33) (RCP9415)

Description	Interleukin 33 (IL-33) is a proinflammatory cytokine belonging to the IL-1 family. IL-33 is cleaved ...
Preservatives	10 mM Tris-HCl pH 8, 50% Glycerol
Form/Appearance	Liquid
Storage	at -20 °C
Tag	N terminal containing a HIS tag
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	22.49 kDa
Uniprot ID	O95760
Expiration Date	6 months from date of receipt.

