

## Product Datasheet

# Radox - Pencillin Binding Protein 2X (PBP-2X) (Western Blot Control) (RCP9490)



<b>Description</b>	Penicillin-binding proteins (PBPs) function in the late steps of murein biosynthesis. Penicillin Bin...
<b>Preservatives</b>	1X Laemmli buffer
<b>Form/Appearance</b>	Liquid. Supplied in 1x Laemmli Buffer (25mM Tris-HCl pH 6.8, 50mM DTT, 1% (w/v) SDS, 0.1% (w/v) Bromophenol Blue, 2.5% Glycerol).
<b>Storage</b>	at -20 °C
<b>Tag</b>	N terminal containing a HIS tag
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA, Western Blotting, etc.
<b>Purity</b>	> 95 %, by SDS-PAGE under reducing conditions
<b>Source</b>	E.coli -derived human recombinant protein
<b>MW</b>	81.18 kDa
<b>Uniprot ID</b>	<b>P59676</b>
<b>Expiration Date</b>	6 months from date of receipt.