

Product Datasheet

Vazyme - Murine RNase Inhibitor (Glycerol-free) (RL301-01)

Description

Murine RNase Inhibitor (Glycerol-free) is a lyophilizable version of Murine RNase Inhibitor (Biorbyt #R301) with excellent thermal stability and storage stability in the glycerol-free state. Compared to the Human RNase Inhibitor, the Murine RNase Inhibitor does not contain oxidation-sensitive cysteines. The superior antioxidant activity makes Murine RNase Inhibitor (Glycerol-free) suitable for experiments (e.g. qPCR assays) that are sensitive to reducing agents. This product does not contain excipients, please add your own as needed.

Storage

Store at -30 ~ -15°C and transport at ≤0°C.

Note

For research use only

Application notes

Supports lyophilization Strong stability

Expiration Date

12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)