

Product Datasheet

Vazyme - AccurSTART Reverse Transcriptase (RW201-01)

Description	AccurSTART Reverse Transcriptase is a directional modified reverse transcriptase based on HiScript I...
Storage	Store at -30 ~ -15°C and transport at ≤0°C
Note	For research use only
Application notes	High specificity: achieves beautiful 'S-shaped' amplification curves. Strong stability A good choice for self-build one-step systems: performs high reverse transcription efficiency in RT-PCR/RT-qPCR systems.
Expiration Date	12 months from date of receipt.

