

Product Datasheet

Vazyme - AccurSTART Reverse Transcriptase (Glycerol-free) (RWL201-01)

| | |
|--------------------------|---|
| Description | AccurSTART Reverse Transcriptase (Glycerol-free) is a lyophilizable version of AccurSTART Reverse Tr... |
| Storage | Store at -30 ~ -15°C and transport at ≤0°C |
| Note | For research use only |
| Application notes | High specificity: achieves beautiful 'S-shaped' amplification curves. Strong stability A good choice for self-build one-step systems: performs high reverse transcription efficiency in RT-PCR/RT-qPCR systems. |
| Expiration Date | 12 months from date of receipt. |

