

## Product Datasheet

# Vazyme - AccurSTART Reverse Transcriptase (Glycerol-free) (RWL201-02)

<b>Description</b>	AccurSTART Reverse Transcriptase (Glycerol-free) is a lyophilizable version of AccurSTART Reverse Tr...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C
<b>Note</b>	For research use only
<b>Application notes</b>	High specificity: achieves beautiful 'S-shaped' amplification curves. Strong stability A good choice for self-build one-step systems: performs high reverse transcription efficiency in RT-PCR/RT-qPCR systems.
<b>Expiration Date</b>	12 months from date of receipt.

