

Product Datasheet

Vazyme - Lipomaster 2000 Transfection Reagent (TL201-02)

Description	Lipomaster 2000 Transfection Reagent is a liposome-mediated transfection reagent. It is applicable f...
Storage	Store a 2 ~ 8°C. Adjust the shipping method according to the destination.
Note	For research use only
Application notes	Broad cell compatibility: Suitable for transfection of various adherent/suspension cells or co-transfection of siRNA and plasmid DNA High transfection efficiency: Efficient transfection for both plasmid DNA and RNA Low cytotoxicity: No need to change medium after transfection
Expiration Date	12 months from date of receipt.

