

Product Datasheet

Vazyme - Lipomaster 3000 Transfection Reagent (TL301-01)

Description Lipomaster 3000 Transfection Reagent is a liposome-mediated transfection reagent. It is applicable f...

Storage Store a 2 ~ 8°C. Adjust the shipping method according to the destination.

Note For research use only

Application notes Superior performance: Efficient transfection of large plasmids and difficult-to-transfect cells Broad compatibility: Only one reagent can be used for transfection of plasmid DNA, mRNA, siRNA, and co-transfection, etc. Keep cell viability: Gentle on cells, low toxicity

Expiration Date 12 months from date of receipt.

