

Product Datasheet

Rat Sult1a1 protein (orb10081)

Description	Recombinant of rat Sult1a1 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10XHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag info: N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged Expression Region: 1-291aa Sequence Info: Full Length Glycerol content: 0.5
Protein Sequence	MEFSRPPLVHVKGIPLIKYFAETIGPLQNFTAWPDDLLISTYPKSGTTWMSEILDMI YQGGKLEKCGRAPIYARVPFLEFKCPGVPSGLETLEETAPAPRLKTHLPLSLLPQSL LDQKVKVYIARNADVVVSYYNFYNMAKLHPDPGTWDSFLENFMDGEVSYGS WYQHVKEWWELRHTHPVLYLFYEDIKENPKREIKKILEFLGRSLPEETVDSIVHHT SFKKMKENCMTNYTTIPTTEIMDHNVSFMRKGTGDKNTFTVAQNERFDAHY AKTMTDCDFKFRCEL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	53.9 kDa
Uniprot ID	P17988
Expiration Date	6 months from date of receipt.