

Product Datasheet

PRRSV-GP5 antibody (orb100850)

Description	PRRSV-GP5 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from PRRSV GP5 (101-200/200 aa)
Target	PRRSV-GP5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q9WJB2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.