
Product Datasheet

Anti-MESP1/APC (orb1009226)

Description

MESP1, also known as bHLHc5, is a 268 amino

helix (b...

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

APC

Tested Applications

ICC, IF

Immunogen

KLH conjugated synthetic peptide derived from human MESP1

Target

MESP1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C

Note

For research use only

Application notes

Optimal dilutions/concentrations should be determined by the end user

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A

MW

30kDa

Uniprot ID

[Q9BRJ9](#)

Dilution Range

ICC 1:50-200, IF 1:50-200

Expiration Date

12 months from date of receipt.