

Product Datasheet

NBS1 p95 (phospho-Ser343) antibody (orb101094)

Description	NBS1 p95 (phospho-Ser343) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NBS1 around the phosphorylation site of Ser343 (SL(p-S)QG)
Target	phospho-NBN (Ser343)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85 kDa
Uniprot ID	O60934
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.