



Product Datasheet

GGT1 antibody (orb101165)


Description	GGT1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human GGT1 (351-460/569 aa)
Target	GGT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	P19440
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
(mo)
Nephroc...



Flow
cytometric
analysis of
Raji cell us...



Flow
cytometric
analysis of
(mo)
Nephroc...