

---

## Product Datasheet

### BMP11 antibody (orb101175)

**Description**

BMP11 antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P, WB

**Immunogen**

KLH conjugated synthetic peptide derived from human BMP11. Please contact us for the exact immunogen sequence. The peptide is available as orb393924.

**Target**

BMP11

**Form/Appearance**

10 mM PBS, 0.02% sodium azide

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**MW**


45 kDa

**Uniprot ID**
[Q9Z217](#), [O95390](#), [Q9Z1W4](#)
**NCBI**
[005811](#), [005802](#)
**Entrez**
[10220](#)
**Dilution Range**


IHC-P: 1:100-600, WB: 1:200-2000, WB: 1:200-1000

**Expiration Date**


12 months from date of receipt.



IHC-P staining of rat colon tissue using...



Immunohistochemical staining of paraffin...



Western blot analysis of mouse colon (la...