



## Product Datasheet

### RAR gamma antibody (orb101688)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | RAR gamma antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human Retinoic Acid Receptor gamma (151-250/454 aa) |
| <b>Target</b>              | RAR gamma   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                             |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                  |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 50 kDa  |
| <b>Uniprot ID</b>          | <b>P13631</b>   |
| <b>Dilution Range</b>      | WB=1:500-2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Immunohistochemical analysis of formalin...



IHC-P of human rectal carcinoma tissue (...)