



Product Datasheet

NRG1 antibody (orb101689)


Description	NRG1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Neuregulin 1 (65-150/640 aa)
Target	HRG beta 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70 kDa
Uniprot ID	Q02297
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100-500, IF=1:200-800, ELISA=1:5000-10000, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse br...



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat brai...