

Product Datasheet

CAMK4 (phospho-Thr196/200) antibody (orb101719)

Description	CAMK4 (phospho-Thr196/200) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CAMKIV around the phosphorylation site of Thr196 + Thr200. (MK(p-T)VCG(p-T)PG)
Target	phospho-CAMK4 (Thr196 + Thr200)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q16566
Expiration Date	12 months from date of receipt.