



Product Datasheet

BCL2 antibody (orb10175)

| | |
|----------------------------|---|
| Description | BCL2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | FC, ICC, IHC-P |
| Immunogen | KLH conjugated synthetic peptide derived from human Bcl-xL (81-200/233 aa) |
| Target | BCL2L1 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 26 kDa |
| Uniprot ID | Q07817 |
| Dilution Range | IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1ug/Test, WB=1:500-2000, IHC-F=1:100-500, IF=1:100, ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of formalin...



Western blot analysis of Rat brain lysat...