



Product Datasheet

SLC14A1 antibody (orb101768)

Description	SLC14A1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human SLC14A1/RACH1 (151-250/389 aa)
Target	SLC14A1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q13336
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Sample: Kidney (Mouse) Lysate at 40 ug P...



Sample: Liver (Mouse) Lysate at 40 ug Pr...