
Product Datasheet

**CD32B (phospho-Tyr292) antibody (FITC)
(orb102270)**

Description

Rabbit Polyclonal to CD32B (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CD32B around the phosphorylation site of Tyr291.
Target	FCGR2B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31kDa
Uniprot ID	P31994
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.