

Product Datasheet

phospho-C-Myc(Thr358)/RBITC antibody (orb1046449)

Description	phospho-C-Myc(Thr358)/RBITC antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	RBITC
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human C-Myc around the phosphorylation site of Thr358
Target	phospho-C-Myc (Thr358)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.