

Product Datasheet

Phospho-STAT3 (Tyr705) Rabbit Polyclonal Antibody (Cy7) (orb1050850)

| | |
|----------------------------|--|
| Description | Phospho-STAT3 (Tyr705) Rabbit Polyclonal Antibody (Cy7) |
| Species/Host | Rabbit |
| Reactivity | Gallus, Human, Mouse, Rat |
| Conjugation | Cy7 |
| Tested Applications | FC, IF |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human STAT3 around the phosphorylation site of Tyr705 AP(p-Y)LK |
| Target | STAT3 |
| Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Antibody Type | Primary Antibody |
| Dilution Range | Flow-Cyt=1µg /test, IF=1:100-500 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)