



Product Datasheet

alpha Defensin 1 antibody (orb10527)

Description	alpha Defensin 1 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human alpha Defensin 1/3 (66-94/94 aa)
Target	DEFA1/alpha Defensin 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6.7 kDa
Uniprot ID	P59665
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human colon tissue using alpha ...