

Product Datasheet

Cyclin D1 (Phospho-Thr286) antibody (orb105916)

Description Rabbit polyclonal to Cyclin D1 (Phospho-Thr286). The protein encoded by this gene belongs to the hig...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC-P, WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of MCF7 whole cell...