


Product Datasheet

JAK2 (Phospho-Tyr931) antibody (orb105936)

Description	Rabbit polyclonal to JAK2 (Phospho-Tyr931). This gene product is a protein tyrosine kinase involved ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...