


Product Datasheet

Chk2 (Phospho-Thr383) antibody (orb105944)

Description	Rabbit polyclonal to Chk2 (Phospho-Thr383). In response to DNA damage and replication blocks, cell C...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...