

## Product Datasheet

### p53 (Phospho-Ser392) antibody (orb105975)

**Description** Rabbit polyclonal antibody to p53 (Phospho-Ser392)

**Species/Host** Rabbit

**Reactivity** Human

**Conjugation** Unconjugated

**Tested Applications** ICC, IF, IHC-P, WB


**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Clonality** Polyclonal

**Dilution Range** WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500

**Expiration Date** 12 months from date of receipt.



Western blot analysis of 293 whole cell ...