


Product Datasheet

HSP27 (Phospho-Ser78) antibody (orb105981)

Description	Rabbit polyclonal to HSP27 (Phospho-Ser78). HSP27 is a small heat shock protein that is regulated bo...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...