

Product Datasheet

PDGFR beta (Phospho-Tyr1021) antibody (orb106022)

Description	Rabbit polyclonal to PDGFR beta (Phospho-Tyr1021). This gene encodes a cell surface tyrosine kinase ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.