

Product Datasheet

Androgen Receptor (Phospho-Ser650) antibody (orb106026)

Description	Rabbit polyclonal to Androgen Receptor (Phospho-Ser650). The androgen receptor gene is more than 90 ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.