

Product Datasheet

B-RAF (Phospho-Thr599) antibody (orb106054)

Description	Rabbit polyclonal to B-RAF (Phospho-Thr599). B-Raf is a tyrosine kinase-like kinase of the RAF famil...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.