

Product Datasheet

PLCG2 (phospho-Tyr753) antibody (orb106071)

| | |
|----------------------------|--|
| Description | PLCG2 (phospho-Tyr753) antibody |
| Species/Host | Rabbit |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IHC-Fr, IHC-P, WB |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 2 around the phosphorylation site of Tyr753 (SL(p-Y)DV) |
| Target | Phospho-PLCG 2 (Tyr753) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 147 kDa |
| Uniprot ID | P16885 |
| Expiration Date | 12 months from date of receipt. |