


Product Datasheet

IGFBP-3 (Phospho-Ser183) antibody (orb106080)

Description	Rabbit polyclonal to IGFBP-3 (Phospho-Ser183). This gene is a member of the insulin-like growth fact...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...