


Product Datasheet

Epo-R (Phospho-Tyr368) antibody (orb106086)

Description	Rabbit polyclonal to Epo-R (Phospho-Tyr368). Erythropoiesis is regulated through the interaction of ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...