


Product Datasheet

NF-kappaB p105/p50 (Phospho-Ser907) antibody (orb106094)

Description	Rabbit polyclonal to NF-kappaB p105/p50 (Phospho-Ser907). NFkB-p105 a transcription factor of the nu...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...