


Product Datasheet

IL-8R beta/CDw128 beta (Phospho-Ser347) antibody (orb106102)

Description	Rabbit polyclonal to IL-8R beta/CDw128 beta (Phospho-Ser347). IL-8R B Receptor for interleukin-8 whi...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 whole c...