

## Product Datasheet

### CDC25A (Phospho-Thr507) antibody (orb106116)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to CDC25A (Phospho-Thr507). CDC25A is a member of the CDC25 family of phosphatases...                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | WB                                                                                                                                      |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>Dilution Range</b>      | WB:1:500-1:2000                                                                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |