


Product Datasheet

CDC25C (Phospho-Ser216) antibody (orb106119)

Description	Rabbit polyclonal to CDC25C (Phospho-Ser216). a member of the MPI phosphatase family. Activates the ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HuvEc whole cel...