


## Product Datasheet

### VEGFR2 (Phospho-Tyr1214) antibody (orb106136)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to VEGFR2 (Phospho-Tyr1214)                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                 |
| <b>Dilution Range</b>      | WB:1:500-1:2000, IHC-P:1:50-1:200                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |



Western blot analysis of HepG2 whole cel...