


Product Datasheet

HNF4 (phospho-Ser313) antibody (orb106148)

Description	HNF4 (phospho-Ser313) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HNF4 alpha around the phosphorylation site of Ser313 (LR(p-S)QV)
Target	phospho-HNF4 (Ser313)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	P41235
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=2ug/test
Expiration Date	12 months from date of receipt.



Flow cytometry of HepG2(black line).Prim...