


Product Datasheet

ZAP-70 (Phospho-Tyr319) antibody (orb106157)

Description	Rabbit polyclonal to ZAP-70 (Phospho-Tyr319). ZAP70 (human) a tyrosine kinase of the Syk family. Ass...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...