


Product Datasheet

JNk1/2/3 (Phospho-Thr183+Tyr185) antibody (orb106161)

Description	Rabbit polyclonal to JNk1/2/3 (Phospho-Thr183+Tyr185). JNK3 a protein kinase of the MAPK family that...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...