

Product Datasheet

SEK1/Mkk4 (Phospho-Thr261) antibody (orb106165)

Description Rabbit polyclonal to SEK1/Mkk4 (Phospho-Thr261). MKK4 dual specificity kinase of the STE7 family tha...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC-P, WB


Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of 293 whole cell ...