

## Product Datasheet

# Tuberin/TSC2 (Phospho-Thr1462) antibody (orb106171)

<b>Description</b>	Rabbit polyclonal to Tuberin/TSC2 (Phospho-Thr1462). Mutations in this gene lead to tuberous scleros...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IHC-P:1:50-1:200
<b>Expiration Date</b>	12 months from date of receipt.